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**EN** **Instructions for Use**  
Sterile, synthetic, non-absorbable surgical suture material

## 1. DESCRIPTION

SABANA PTFE suture is monofilament, synthetic, non-absorbable surgical suture composed of 100% polytetrafluoroethylene without any additives. The molecular formula is (C<sub>2</sub>F<sub>4</sub>)<sub>n</sub>. SABANA PTFE suture is undyed. SABANA PTFE sutures are available in a range of gauge sizes and lengths, non-needle or attached to stainless steel needles of varying types and sizes. SABANA PTFE sutures comply with the requirements of the European Pharmacopoeia for Sterile Non-Absorbable Strands and the requirements of the United States Pharmacopoeia for Non-Absorbable Surgical Sutures.

## 2. INDICATIONS

SABANA PTFE suture is indicated for use in all types of soft tissue approximation and/or ligation, including cardiovascular, dental and general surgeries, as well as repair of the dura mater. The device is not indicated for use in ophthalmic surgery, microsurgery and peripheral neural tissue.

## Application

SABANA PTFE sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

## Material Performance

SABANA PTFE sutures are not absorbed. Due to the relative biological inertness of PTFE sutures, they are recommended for use where the least possible suture reaction is desired.

## 3. CONTRAINDICATIONS

None known.

## 4. WARNINGS / PRECAUTIONS / INTERACTIONS

- Do not use if package is opened or damaged. Discard opened unused sutures.
- Do not use after expiry date.
- The safety and effectiveness of this suture in ophthalmic, microsurgical, and peripheral neural applications has not been established.
- Tissue invasion of SABANA PTFE suture can result in attachment of the suture to the tissue it penetrates in long-term use. Such attachment may make removal of the suture difficult.
- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing SABANA PTFE sutures for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.
- Acceptable surgical practice should be followed for the management of infected or contaminated wounds.
- As with any exogenous material, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation.
- Adequate knot security requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting SABANA PTFE sutures.
- When handling this or any other suture, care should be taken to avoid damage. Avoid crushing or crimping damage due to the application of surgical instruments such as forceps or needle holders.
- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury.
- Discard used needles in "Sharps" container.
- Dispose of material in accordance with all the state, local, and hospital regulations. Responsibility for proper waste disposal is with the owner of the waste.
- Do not reuse: Infection hazard for patients and/or users and impairment of products functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.
- Do not resterilize: Infection hazard for patients and/or users and impairment of products functionality due to use of resterilized suture. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

## 5. ADVERSE REACTIONS

Adverse reactions associated with the use of this device include minimal initial inflammatory tissue reaction and transient local irritation at the wound site. As with all foreign bodies, SABANA PTFE sutures may exacerbate an existing infection.

## 6. STERILITY

SABANA PTFE sutures are sterilized by ethylene oxide gas. Do not resterilize!  
Do not use if package is opened or damaged! Discard opened unused sutures!

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## 7. STORAGE CONDITIONS

Recommended storage conditions for SABANA PTFE: store at a temperature below 25°C, away from moisture, corrosion and direct heat. Do not use after expiry date.

## 8. LABELING SYMBOLS

 yyyy-mm Date of Manufacture (year, month)

 Undyed, monofilament, non-absorbable suture material

 LOT Batch designation

 STERILE EO Sterilized with ethylene oxide

 Comply with the instructions for use.

 Caution, consult accompanying document

 Do not reuse.

 Do not resterilize.

 yyyy-mm Use by date (year-month)

 Store in a dark location away from sunlight.

 Keep dry

 25°C Upper temperature limit for storage

 Do not use if the packaging is damaged.

**CE** 0123 CE mark with code number of the Notified Body

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