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DESCRIPTION

SABAsorb® sutures are synthetic, absorbable, sterile surgical sutures composed of Polyglycolic acid (PGA). The empirical formula of the polymer is (C₃H₄O₂)_n. SABAsorb® sutures are available undyed and dyed violet with D&C Violet No. 2 (Colour Index number 60725).

SABAsorb® sutures are available as braided strands in EP sizes 0.4-6 (USP sizes 8-0 through 3 or 4), in a variety of lengths, with and without stainless steel needles of varying types and sizes. Braided sutures are uniformly coated with polycaprolactone and calcium stearate.

SABAsorb® sutures are also available as monofilament in EP sizes 0.2-0.3 (USP size 10-0 through 9-0), in a variety of lengths, with and without stainless steel needles of varying types and sizes. Monofilament sutures are uncoated.

SABAsorb® suture complies with the requirements of the European Pharmacopoeia for "Sutures, Sterile Synthetic Absorbable Braided" and the requirements of United States Pharmacopoeia for "Absorbable Surgical Suture".

INDICATIONS

SABAsorb® sutures are intended for use in general soft tissue approximation and/or ligation.

Application

SABAsorb® sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

Performance

SABAsorb® sutures elicits a minimal initial inflammatory reaction in tissues and are eventually replaced with an ingrowth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of sutures occurs by means of hydrolysis, where the polymer degrades to glycolic acid which is subsequently absorbed and eliminated by the body. Absorption begins as a loss of tensile strength followed by a loss of mass. Implantation studies in rats show the following profile.

Days implantation	Approximate % original strength remaining
14 days	75 %
21 days	40 %

Absorption of SABAsorb® sutures is essentially complete between 60 and 90 days.

CONTRAINDICATIONS

These sutures, being absorbable should not be used where extended approximation of tissues under stress is required. SABAsorb® sutures are not intended for use in cardiovascular tissue, neurosurgery or ophthalmic surgery.

WARNINGS/PRECAUTIONS/INTERACTIONS

- Do not use if package is opened or damaged. Discard opened unused sutures.
- Do not use after expiry date.
- Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing SABAsorb® suture for wound closure, as risk of wound dehiscence, may vary with the site of application and the suture material used. Surgeons should consider the in vivo performance (described in PERFORMANCE section) when selecting a suture.
- As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tract may result in calculus formation. As an absorbable suture, it may act transiently as a foreign body.
- Acceptable surgical practice should be followed for the management of contaminated or infected wounds.
- As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distension, or which may require additional support.
- Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated.
- Under some circumstances, notably orthopaedic procedures, immobilisation of joints by external support may be employed at the discretion of the surgeon.
- Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.
- Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.
- This suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.
- When handling this or any other suture, care should be taken to avoid damage. Avoid crushing or crimping damage due to the application of surgical instruments such as forceps or needle holders.
- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Broken needles may result in extended or additional surgeries or residual foreign bodies.
- Adequate knot security requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting any monofilament suture.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens.

- Discard used needles in "Sharps" container.
- Dispose of material in accordance with all the state, local, and hospital regulations. Responsibility for proper waste disposal is with the owner of the waste.
- Do not reuse: Infection hazard for patients and/or users and impairment of products functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.
- Do not resterilize: Infection hazard for patients and/or users and impairment of products functionality due to use of resterilized suture. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

ADVERSE REACTIONS

Adverse effects associated with the use of device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens.

STERILITY

SABAsorb® sutures are sterilized by ethylene oxide gas. Do not resterilize! Do not use if package is opened or damaged. Discard opened unused sutures.

STORAGE

Recommended storage conditions: stored at a temperature below 25°C, away from moisture corrosion and direct heat.

SYMBOLS USED ON LABELLING

Sterilized using ethylene oxide

Do not resterilize

Do not reuse

Do not use if package is damaged

Caution, consult accompanying document

Consult instructions for use

Keep away from sunlight

Temperature limitation

Keep dry

Catalogue number

Batch code

Date of manufacture

Used by date

CE-Mark and Identification number of Notified Body

Manufacturer

Authorised representative in the European community

Fragile, handle with care

IMPORTER AND DISTRIBUTOR

Sabana Medizinbedarf GmbH
Alwinenstr. 2
65189 Wiesbaden / Germany
Phone: +49 (0) 6 11 - 50 46 40 8-0
Fax: +49 (0) 6 11 - 50 46 40 8-9
E-Mail: info@sabana.de
www.sabana.de

MANUFACTURER

Foosin Medical Supplies Inc., Ltd.
No. 20, Xingshan Road, Weihai Torch Hi-tech Science Park
264210 Weihai, Shandong Province
People's Republic of China
Phone: +86-631-5661991
Fax: +86-631-5661996

EC REPRESENTATIVE

MedNet EC-REP GmbH
Borkstr. 10
48163 Muenster / Germany