DE

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HERSTELLER

Foosin Medical Supplies Inc., Ltd.
No. 20, Xingshan Road, Weihai Torch Hi-tech Science Park
264210 Weihai, Shandong Province
People's Republic of China
Telefon: +86 – 6 31 – 5 66 19 91
Telefax: +86 – 6 31 – 5 66 19 96

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EN

Instructions for Use

Sterile, synthetic absorbable sutures (Polyglycolic acid)

DESCRIPTION

SABAsorb rapid® sutures are synthetic, absorbable, braided, sterile surgical sutures composed of Polyglycolic acid (PGA). The empirical formula of the polymer is (C₃H₂O₂)_n. The characteristic rapid loss of strength is achieved by use of a polymer material with a lower molecular weight than regular SABAsorb® suture. SABAsorb rapid® sutures are available undyed and dyed violet with D&C Violet No. 2 (Colour Index number 60725). SABAsorb rapid® sutures are coated with polycaprolactone and calcium stearate. SABAsorb rapid® sutures are available in a range of gauge sizes and lengths, with and without stainless steel needles of varying types and sizes. SABAsorb rapid® sutures comply with the requirements of the European Pharmacopoeia for “Sutures, Sterile Synthetic Absorbable Braided” and the requirements of United States Pharmacopoeia for “Absorbable Surgical Suture” (except for an occasional slight oversize in diameter).

EP Gauge	Maximum Oversize (mm)
0.7	0.026
1.0	0.026
1.5	0.026
2	0.026
3	0.026
3.5	0.051
4	0.051
5	0.051

INDICATIONS

SABAsorb rapid® sutures are intended for use in general soft tissue approximation where only short term wound support is required, including ophthalmic (e.g. conjunctiva) procedures.

Application

SABAsorb rapid® sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

Performance

SABAsorb rapid® sutures elicits a minimal initial inflammatory reaction in tissues and are eventually replaced with an in-growth of fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of sutures occurs by means of hydrolysis, where the polymer degrades to glycolic acid which is subsequently absorbed and eliminated by the body.

Absorption begins as a loss of tensile strength followed by a loss of mass. Implantation studies in rats show the following profile.

Days implantation	Approximate % original strength remaining
7 days	55 %
14 days	20 %
21 days	5 %

Absorption of SABAsorb rapid® sutures are essentially complete between 42 and 63 days.

CONTRAINDICATIONS

Due to the rapid loss of tensile strength, SABAsorb rapid® should not be used where extended approximation of tissues under stress is required or where wound support or ligation beyond 7 days is required. SABAsorb rapid® suture is not for use in cardiovascular and neurological tissues.

WARNINGS/PRECAUTIONS/INTERACTIONS

- Do not use if package is opened or damaged. Discard opened unused sutures.
- Do not use after expiry date.
- Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing SABAsorb rapid® for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Surgeons should consider the in vivo performance (described in PERFORMANCE section) when selecting a suture.
- As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tract may result in calculus formation. As an absorbable suture, SABAsorb rapid® may act transiently as a foreign body.
- Acceptable surgical practice should be followed for the management of contaminated or infected wounds.
- As this is an absorbable suture material, the use of supplemental non-absorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distension, or which may require additional support.
- Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated.
- Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.
- Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.
- This suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.
- When handling this or any other suture material, care should be taken to avoid damage. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (⅓) to one-half (½) of the distance from the attachment end to the point. Grasping in the point area could impair

the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Broken needles may result in extended or additional surgeries or residual foreign bodies.

- Adequate knot security requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstances and the experience of the surgeon.
- Users should exercise caution when handing surgical needles to avoid inadvertent needle stick injury. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens.
- Discard used needles in “Sharps” container.
- Dispose of material in accordance with all the state, local, and hospital regulations. Responsibility for proper waste disposal is with the owner of the waste.
- Do not reuse: Infection hazard for patients and/or users and impairment of products functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.
- Do not resterilize: Infection hazard for patients and/or users and impairment of products functionality due to use of resterilized suture. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

ADVERSE REACTIONS

Adverse effects associated with the use of device include wound dehiscence, failure to provide adequate wound support in closure of the sites where expansion, stretching or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and transitory local irritation at the wound site.

Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens.

STERILITY

SABAsorb rapid® sutures are sterilized by ethylene oxide gas. Do not resterilize! Do not use if package is opened or damaged. Discard opened unused sutures!

STORAGE

Recommended storage conditions: stored at a temperature below 25°C, away from moisture corrosion and direct heat.

SYMBOLS USED ON LABELLING

 Sterilized using ethylene oxide

 Do not resterilize

 Do not reuse

 Do not use if package is damaged

 Caution, consult accompanying document

 Consult instructions for use

 Keep away from sunlight

 Temperature limitation

 Keep dry

 Catalogue number

 Batch code

 Date of manufacture

 Used by date

 CE-Mark and Identification number of Notified Body

 Manufacturer

 Authorised representative in the European community

 Fragile, handle with care

IMPORTER AND DISTRIBUTOR

Sabana Medizinbedarf GmbH
Alwinenstr. 2
65189 Wiesbaden / Germany
Phone: +49 (0) 6 11 – 50 46 40 8-0
Fax: +49 (0) 6 11 – 50 46 40 8-9
E-Mail: info@sabana.de
www.sabana.de

MANUFACTURER

Foosin Medical Supplies Inc., Ltd.
No. 20, Xingshan Road, Weihai Torch Hi-tech Science Park
264210 Weihai, Shandong Province
People's Republic of China
Phone: +86 – 6 31 – 5 66 19 91
Fax: +86 – 6 31 – 5 66 19 96

EC REPRESENTATIVE

MedNet EC-REP GmbH
Borkstr. 10
48163 Muenster / Germany